

## Product datasheet for **RC205673L4V**

### IGF2BP2 (NM\_006548) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	IGF2BP2 (NM_006548) Human Tagged ORF Clone Lentiviral Particle
Symbol:	IGF2BP2
Synonyms:	IMP-2; IMP2; VICKZ2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006548
ORF Size:	1794 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205673).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006548.4</a>
RefSeq Size:	3676 bp
RefSeq ORF:	1800 bp
Locus ID:	10644
UniProt ID:	<a href="#">Q9Y6M1</a>
Cytogenetics:	3q27.2
Domains:	RRM, KH
MW:	66 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a protein that binds the 5' UTR of insulin-like growth factor 2 (IGF2) mRNA and regulates its translation. It plays an important role in metabolism and variation in this gene is associated with susceptibility to diabetes. Alternative splicing and promoter usage results in multiple transcript variants. Related pseudogenes are found on several chromosomes. [provided by RefSeq, Sep 2016]