

Product datasheet for **RC205682L3V**

MSX1 (NM_002448) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MSX1 (NM_002448) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MSX1
Synonyms:	ECTD3; HOX7; HYD1; STHAG1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_002448
ORF Size:	909 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205682).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002448.3 , NP_002439.2
RefSeq Size:	1940 bp
RefSeq ORF:	912 bp
Locus ID:	4487
UniProt ID:	P28360
Cytogenetics:	4p16.2
Domains:	homeobox
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

MW: 31.5 kDa

Gene Summary: This gene encodes a member of the muscle segment homeobox gene family. The encoded protein functions as a transcriptional repressor during embryogenesis through interactions with components of the core transcription complex and other homeoproteins. It may also have roles in limb-pattern formation, craniofacial development, particularly odontogenesis, and tumor growth inhibition. Mutations in this gene, which was once known as homeobox 7, have been associated with nonsyndromic cleft lip with or without cleft palate 5, Witkop syndrome, Wolf-Hirschorn syndrome, and autosomal dominant hypodontia. [provided by RefSeq, Jul 2008]