

Product datasheet for **RC205690L1V**

PLAGL2 (NM_002657) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PLAGL2 (NM_002657) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PLAGL2
Synonyms:	ZNF900
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_002657
ORF Size:	1488 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205690).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002657.2
RefSeq Size:	5640 bp
RefSeq ORF:	1491 bp
Locus ID:	5326
UniProt ID:	Q9UPG8
Cytogenetics:	20q11.21
Domains:	zf-C2H2
Protein Families:	Transcription Factors



[View online »](#)

MW: 54.4 kDa

Gene Summary: Pleiomorphic adenoma gene-like 2 is a zinc-finger protein that recognizes DNA and/or RNA. [provided by RefSeq, Jul 2008]