

Product datasheet for RC205726L1

VKORC1L1 (NM_173517) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	VKORC1L1 (NM_173517) Human Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	VKORC1L1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205726).
Restriction Sites:	Sgfl-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_173517
ORF Size:	528 bp




[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_173517.3
RefSeq Size:	5913 bp
RefSeq ORF:	531 bp
Locus ID:	154807
UniProt ID:	Q8N0U8
Cytogenetics:	7q11.21
Protein Families:	Transmembrane
MW:	19.8 kDa
Gene Summary:	This gene encodes an enzyme important in the vitamin K cycle, which is involved in the carboxylation of glutamate residues present in vitamin K-dependent proteins. The encoded enzyme catalyzes the de-epoxidation of vitamin K 2,3-epoxide. Oxidative stress may upregulate expression of this gene and the encoded protein may protect cells and membrane proteins from oxidative damage. This gene and a related gene (Gene ID: 79001) may have arisen by gene duplication of an ancestral gene. [provided by RefSeq, Oct 2016]

Product images:



Circular map for RC205726L1



Double digestion of RC205726L1 using SgfI and MluI