

Product datasheet for RC205732A1V

Human C2orf15 (NM_144706) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human C2orf15 (NM_144706) AAV Particle
Tag: Myc-DDK
Symbol: C2orf15
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC205732 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGTCCTCAACTCTGGGGAAGTTAAGTAATCAAGTTGAAGAAACACTTCCACTACTTAAAAAGGTACCTG
 CAAATTACTTTACATTTGTTTCAGCTATCCTAATGGGATTTTCACTTAGTAAATCTGCTACTCAGGTATC
 TGCTATACATATGGATTCAAAAGTGGATGATCACTTAATACGAGGGACTGAAAAAGCAGTTGGAACCA
 GCGACTCAGTTATTTCAAAACACCAAGAAAATAAGATTAGAAGACACAAATCAAGAAAACCTTTACAAGGA
 TTGAAGGGACTGGCAGGATCTCTTTCTGGGAAAGCCTTGGGTTCAAGTGTATATGTCAAAGAAAGTGA
 TGGACTAGAAATGACAGATGTGGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC205732 protein sequence
 Red=Cloning site Green=Tags(s)

MSSTLGKLSNQVEETLPLKKVPANYFHICSAILMGFSLSKSATQVSAIHMSKVDHDLIRGTEKSRLEP
 ATQLFQNTKKIRLEDTNQENFTRIEGTGTGSLSGKALGSVVYVKESDGLMTDVE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_144706
ORF Size: 375 bp


[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_144706.2</u>
RefSeq Size:	1430 bp
RefSeq ORF:	276 bp
Locus ID:	150590
UniProt ID:	<u>Q8WU43</u>
Cytogenetics:	2q11.2
MW:	13.7 kDa