

Product datasheet for **RC205748A1V**

Human CUTC (NM_015960) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CUTC (NM_015960) AAV Particle
Tag:	Myc-DDK
Symbol:	CUTC
Synonyms:	CGI-32
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC205748 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAAAAGGCAGGGGGCCTCCTCTGAGCGAAAACGAGCGCGGATACCGTCCGGGAAGCCGGAGCAGCAA
ATGGATTTCTCATGGAAGTTTGTGTGATTGAGTGAATCAGCTGTGAATGCAGAAAGAGGAGGTGCTGA
TCGGATTGAATTATGTTCTGGTTTATCAGAGGGGGGAACACCCAGCATGGGTGTCTTCAAGTAGTG
AAGCAGAGTGTTTCAGATCCAGTTTTTGTGATGATTCGCCACGGGGAGGTGATTTTTTGTATTCAGATC
GTGAAATTGAGGTGATGAAGGCTGACATTCGTCTTGCCAAGCTTTATGGTGCTGATGGTTTGGTTTTGG
GGCATTGACTGAAGATGGACACATTGACAAAGAGCTGTGTATGTCCTTATGGCTATTTGCCGCCCTCTG
CCAGTCACTTTCCACCGAGCCTTTGACATGGTTCATGATCCAATGGCAGCTCTGGAGACCCTCTTAACCT
TGGGATTTGAACGCGTGTGACCAAGTGGATGTGACAGTTCAGCATTAGAAGGGCTACCCCTAATAAAGCG
ACTCATTGAGCAGGCAAAGGCAGGATTGTGTAATGCCAGGAGGTGGTATAACAGACAGAAATCTACAA
AGGATCCTTGAGGGTTCAGGTGCTACAGAATCCACTGTTCTGCTCGGTCTACTAGAGACTCGGGAATGA
AGTTTCGAAATTCATCTGTTGCCATGGGAGCCTCACTTTCTTGCTCAGAATATCCCTAAAGGTAACAGA
TGTGACCAAAGTAAGGACTTTGAATGCTATCGCAAAGAATCCTCGTG

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC205748 protein sequence Red=Cloning site Green=Tags(s) MKRQGASSERKRARIPSGKAGAANGFLMEVCVDSVESAVNAERGGADRIELCSGLSEGGTTPSMGVLQVV KQSVQIPVFVMIRPRGGDFLYSDREIEVMKADIRLAKLYGADGLVFGALTEDGHIDKELCMSLMAICRPL PVTFFHRAFDMVHDPMAALETLLTLGFERVLTSGCDSSALEGLPLIKRLIEQAKGRIVVMPGGGITDRNLQ RILEGSGATEFHCSARSTRDSGMKFRNSSVAMGASLSCSEYSLKVTDVTKVRTLNIAIAKNILV SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_015960
ORF Size:	819 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_015960.1</u>
RefSeq Size:	1378 bp
RefSeq ORF:	822 bp
Locus ID:	51076
UniProt ID:	<u>Q9NTM9</u>
Cytogenetics:	10q24.2
MW:	29.3 kDa