

## Product datasheet for **RC205748L4V**

### CUTC (NM\_015960) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CUTC (NM_015960) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CUTC
Synonyms:	CGI-32
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_015960
ORF Size:	819 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205748).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_015960.1</a>
RefSeq Size:	1378 bp
RefSeq ORF:	822 bp
Locus ID:	51076
UniProt ID:	<a href="#">Q9NTM9</a>
Cytogenetics:	10q24.2
MW:	29.3 kDa



[View online »](#)

**Gene Summary:**

Members of the CUT family of copper transporters are associated with copper homeostasis and are involved in the uptake, storage, delivery, and efflux of copper (Gupta et al., 1995 [PubMed 7635807]; Li et al., 2005 [PubMed 16182249]).[supplied by OMIM, Mar 2008]