

Product datasheet for **RC205751L3V**

FXR1 (NM_001013438) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FXR1 (NM_001013438) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FXR1
Synonyms:	FXR1P; MYOPMIL; MYORIBF
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001013438
ORF Size:	1617 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205751).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001013438.1
RefSeq Size:	8475 bp
RefSeq ORF:	1620 bp
Locus ID:	8087
UniProt ID:	P51114
Cytogenetics:	3q26.33
Protein Families:	Druggable Genome
MW:	60.8 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene is an RNA binding protein that interacts with the functionally-similar proteins FMR1 and FXR2. These proteins shuttle between the nucleus and cytoplasm and associate with polyribosomes, predominantly with the 60S ribosomal subunit. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]