

Product datasheet for **RC205764A1V**

Human FBLIM1 (NM_017556) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human FBLIM1 (NM_017556) AAV Particle
Tag:	Myc-DDK
Symbol:	FBLIM1
Synonyms:	CAL; FBLP-1; FBLP1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC205764 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCTCAAAGCCTGAGAAGAGGGTGGCATCGTCTGTCTTTATCACCCTGGCACCCCCGCGCCGCGATG
TGGCCGTGGCGGAGGAAGTGAGGCAGGCAGTTTGTGAGGCCGGCGTGGCCGCCCTGGGAGGCTCCTGC
CCCCATGAAGACACCCGAGGCTGGCTTGGCGGGAGGCCAGCCCTGGACAACCCCTGGCAGAGCTGCA
GCCACAGTGCCGGCTGCACCTATGCAGCTCTTCAATGGAGGATGCCACCCCTCCTCCTGTCTGGATG
GTGAGGACGTGCTTCTGACCTGGACCTCTCCACCCCTCCACCGCCCTCCAGTGCTTCTGCCTTC
TGAAGAGGAGGCTCCTGCTCCAATGGGGCCTCACTCATTGCAGACTTAGAGCAGCTGCACCTGTCCCG
CCCCGCCCCCACCACAGGCCCCAGCGAGGGACCTTCAGTCCAGCCCGTCCCTCAGGCCCATGGAGG
AAGAGCTGCCACCTCCCCCGGCAGAACCTGTTGAGAAAGGGGCATCCACAGACATCTGTGCCTTCTGCCA
CAAGACCGTGTCCCCGAGAGCTGGCTGTGGAGGCCATGAAGAGGCAGTACCATGCCAGTGCTTCACG
TGCCGCACCTGCCGCCGACGCTGGCTGGGAGAGCTTCTACCAGAAGGATGGGCGACCCCTCTGCGAAC
CCTGCTACCAGGACACACTGGAGAGGTGCGCAAGTGTGGCGAGGTGGTCCGGGACCACATCATAGGGC
CCTGGGCCAGGCTTCCACCCCTCCTGCTTACGTGTGTGACCTGCGCCCGGTGCATTGGGGATGAGAGC
TTTGCCCTGGGCAGCCAGAACGAGGTGTACTGCCTGGACGACTTCTACAGGAAATTCGCCCCGCTGCA
GCATCTGTGAAAATCCCATCATCCCTCGGGATGGGAAAGATGCCTTCAAATCGAATGCATGGGAAGAAA
CTTCCATGAAAATTGCTACAGGTGTGAGGACTGCAGGATCCTCCTGTCTGTGCGAGCCACGGACCAAGGC
TGCTACCCCTGAACAACCATCTCTTCTGCAAGCCATGCCATGTGAAGCGGAGTGCTGCGGGGTGCTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



View online »

Protein Sequence:	>RC205764 protein sequence Red=Cloning site Green=Tags(s) MASKPEKRVASSVFITLAPPRRDVAVAEVVRQAVCEARRGRPWEAPAPMKTPEAGLAGRPSWTTTPGRAA ATVPAAPMQLFNGGCPPPPVLGDGDLPLDLLPPPPPPPVLLPSEEEAPAPMGASLIADLEQLHLSP PPPPQAPAEGPSVQGPLRPMEEELPPPPAEPVEKGASTDICAFCCHKTVFPRELAVEAMKRQYHAQCFT CRTCRRQLAGQSFYQKDGRLCEPCYQDTLERCGKCGEVVRDHIIRALGQAFHPSCFTCVTCARCIGDES FALGSQNEVYCLDDFYRKFAVCSICENPIIPRDGKDAFKIECMGRNFHENCYRCEDCRILLSVEPTDQG CYPLNNHLFCKPCHVKRSAAGCC TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_017556
ORF Size:	1119 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_017556.2
RefSeq Size:	3363 bp
RefSeq ORF:	1122 bp
Locus ID:	54751
UniProt ID:	Q8WUP2
Cytogenetics:	1p36.21
MW:	40.7 kDa