

## Product datasheet for **RC205765L1V**

### **SNX21 (NM\_152897) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	SNX21 (NM_152897) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SNX21
Synonyms:	C20orf161; dj337O18.4; PP3993; SNX-L; SNXL
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_152897
ORF Size:	597 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205765).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_152897.1</a>
RefSeq Size:	1921 bp
RefSeq ORF:	600 bp
Locus ID:	90203
UniProt ID:	<a href="#">Q969T3</a>
Cytogenetics:	20q13.12
Protein Families:	Druggable Genome
MW:	21.9 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein does not contain a coiled coil region, like some family members. The specific function of this protein has not been determined. Multiple transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]