

## Product datasheet for **RC205787L3V**

### SNIP1 (NM\_024700) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SNIP1 (NM_024700) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SNIP1
Synonyms:	PML1; PMRED
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_024700
ORF Size:	1188 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205787).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024700.2</a> , <a href="#">NP_078976.2</a>
RefSeq Size:	4605 bp
RefSeq ORF:	1191 bp
Locus ID:	79753
UniProt ID:	<a href="#">Q8TAD8</a>
Cytogenetics:	1p34.3
Domains:	FHA
MW:	45.8 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a protein that contains a coiled-coil motif and C-terminal forkhead-associated (FHA) domain. The encoded protein functions as a transcriptional coactivator that increases c-Myc activity and inhibits transforming growth factor beta (TGF-beta) and nuclear factor kappa-B (NF-kB) signaling. The encoded protein also regulates the stability of cyclin D1 mRNA, and may play a role in cell proliferation and cancer progression. Mutations in this gene are a cause of psychomotor retardation, epilepsy, and craniofacial dysmorphism (PMRED). [provided by RefSeq, Mar 2012]