

Product datasheet for **RC205792L3V**

ABHD12B (NM_181533) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ABHD12B (NM_181533) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ABHD12B
Synonyms:	BEM46L3; C14orf29; c14_5314
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_181533
ORF Size:	765 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205792).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_181533.3 , NP_853511.2
RefSeq Size:	1580 bp
RefSeq ORF:	768 bp
Locus ID:	145447
UniProt ID:	Q7Z5M8
Cytogenetics:	14q22.1
MW:	28.6 kDa



[View online »](#)