

Product datasheet for **RC205795L4V**

DLX3 (NM_005220) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DLX3 (NM_005220) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DLX3
Synonyms:	AI4; TDO
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_005220
ORF Size:	861 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205795).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005220.2
RefSeq Size:	2613 bp
RefSeq ORF:	864 bp
Locus ID:	1747
UniProt ID:	O60479
Cytogenetics:	17q21.33
Domains:	homeobox
Protein Families:	Druggable Genome, Transcription Factors




[View online »](#)

MW: 31.7 kDa

Gene Summary: Many vertebrate homeo box-containing genes have been identified on the basis of their sequence similarity with Drosophila developmental genes. Members of the Dlx gene family contain a homeobox that is related to that of Distal-less (Dll), a gene expressed in the head and limbs of the developing fruit fly. The Distal-less (Dlx) family of genes comprises at least 6 different members, DLX1-DLX6. Trichodontoosseous syndrome (TDO), an autosomal dominant condition, has been correlated with DLX3 gene mutation. This gene is located in a tail-to-tail configuration with another member of the gene family on the long arm of chromosome 17. Mutations in this gene have been associated with the autosomal dominant conditions trichodontoosseous syndrome and amelogenesis imperfecta with taurodontism. [provided by RefSeq, Jul 2008]

Product images:



[RC205795L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC205795L4V particle to overexpress human DLX3-mGFP fusion protein.