

Product datasheet for RC205797L4

https://www.origene.com

techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

PLTP (NM_006227) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: PLTP (NM_006227) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: **PLTP**

BPIFE; HDLCQ9 Synonyms:

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

Sgfl-Mlul

E. coli Selection: Chloramphenicol (34 ug/mL)

The ORF insert of this clone is exactly the same as(RC205797). **ORF Nucleotide**

Sequence:

Restriction Sites: Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_006227

ORF Size: 1479 bp





PLTP (NM_006227) Human Tagged Lenti ORF Clone - RC205797L4

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 006227.2</u>

 RefSeq Size:
 2100 bp

 RefSeq ORF:
 1482 bp

 Locus ID:
 5360

 UniProt ID:
 P55058

Cytogenetics: 20q13.12

Domains: BPI1

Protein Families: Druggable Genome, Secreted Protein

Protein Pathways: PPAR signaling pathway

MW: 54.7 kDa


Gene Summary: The protein encoded by this gene is one of at least two lipid transfer proteins found in human

plasma. The encoded protein transfers phospholipids from triglyceride-rich lipoproteins to high density lipoprotein (HDL). In addition to regulating the size of HDL particles, this protein may be involved in cholesterol metabolism. At least two transcript variants encoding different

isoforms have been found for this gene. [provided by RefSeq, Jul 2008]



Product images:



Circular map for RC205797L4