

Product datasheet for **RC205837A1V**

Human MRPS16 (NM_016065) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human MRPS16 (NM_016065) AAV Particle
Tag: Myc-DDK
Symbol: MRPS16
Synonyms: CGI-132; COXPD2; MRP-S16; RPMS16
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC205837 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGTCCACCTCACTACTCTCTCTGCAAGGCCTACCGTGGGGGCCACTTAACCATCCGCCTTGGCCTGG
GTGGCTGCACCAATCGGCCGTTCTACCGCATTGTGGCTGCTCACACAAGTGTCACAGGGATGGCCGTTT
CGTAGAGCAGCTGGGCTCCTATGATCCATTGCCAACAGTCATGGAGAAAACTCGTTGCCCTCAACCTA
GACAGGATCCGTCATTGGATTGGCTGCGGGGCCACCTCTCTAAGCCTATGGAAAAGCTTCTGGGTCTTG
CTGGCTTTTTCCCTCTGCATCCTATGATGATCACAATGCTGAGAGACTGCGAAGGAAACGGGCACGTGA
AGTCCTGTAGCTTCTCAGAAAACAGATGCAGAAGCTACAGATACAGAGGCTACAGAAACA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC205837 protein sequence
Red=Cloning site Green=Tags(s)

MVHLTLLCKAYRGHLLTIRLALGGCTNRPFYRIVAHNKCPDRGRFVEQLGSYDPLPNSHGEKLVALLNL
DRIRHWIGCGAHLKPMKLLGLAGFFPLHPMMITNAERLRRKRAREVLLASQKTDAEATDTEATET

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_016065



[View online »](#)

ORF Size:	411 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_016065.3</u>
RefSeq Size:	2651 bp
RefSeq ORF:	414 bp
Locus ID:	51021
UniProt ID:	<u>Q9Y3D3</u>
Cytogenetics:	10q22.2
MW:	15.3 kDa