

Product datasheet for **RC205838L3V**

AMDHD2 (NM_015944) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	AMDHD2 (NM_015944) Human Tagged ORF Clone Lentiviral Particle
Symbol:	AMDHD2
Synonyms:	CGI-14
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_015944
ORF Size:	1317 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205838).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015944.2
RefSeq Size:	1604 bp
RefSeq ORF:	1320 bp
Locus ID:	51005
UniProt ID:	Q9Y303
Cytogenetics:	16p13.3
Domains:	Amidohydro_1
Protein Pathways:	Amino sugar and nucleotide sugar metabolism



[View online »](#)

MW: 46.9 kDa

Gene Summary: Hydrolyzes the N-glycolyl group from N-glycolylglucosamine 6-phosphate (GlcNGc-6-P) in the N-glycolylneuraminic acid (Neu5Gc) degradation pathway. Although human is not able to catalyze formation of Neu5Gc due to the inactive CMAHP enzyme, Neu5Gc is present in food and must be degraded.[UniProtKB/Swiss-Prot Function]