

Product datasheet for **RC205848L4V**

MSH2 (NM_000251) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MSH2 (NM_000251) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MSH2
Synonyms:	COCA1; FCC1; hMSH2; HNPCC; HNPCC1; LCFS2; MMRC52
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_000251
ORF Size:	2802 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205848).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000251.1
RefSeq Size:	3226 bp
RefSeq ORF:	2805 bp
Locus ID:	4436
UniProt ID:	P43246
Cytogenetics:	2p21-p16.3
Domains:	MutS_V, MutS_I, MutS_III, MutS_II, MutS_IV
Protein Families:	Druggable Genome, Stem cell - Pluripotency



[View online »](#)

Protein Pathways: Colorectal cancer, Mismatch repair, Pathways in cancer

MW: 104.7 kDa

Gene Summary: This locus is frequently mutated in hereditary nonpolyposis colon cancer (HNPCC). When cloned, it was discovered to be a human homolog of the E. coli mismatch repair gene mutS, consistent with the characteristic alterations in microsatellite sequences (RER+ phenotype) found in HNPCC. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2012]