

## Product datasheet for **RC205852L4V**

### Ribophorin II (RPN2) (NM\_002951) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ribophorin II (RPN2) (NM_002951) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Ribophorin II
Synonyms:	RIBIIR; RPN-II; RPNII; SWP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_002951
ORF Size:	1893 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205852).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_002951.2</a>
RefSeq Size:	2538 bp
RefSeq ORF:	1896 bp
Locus ID:	6185
UniProt ID:	<a href="#">P04844</a>
Cytogenetics:	20q11.23
Protein Families:	Transmembrane
Protein Pathways:	Metabolic pathways, N-Glycan biosynthesis



[View online »](#)

**MW:** 69.3 kDa

**Gene Summary:** This gene encodes a type I integral membrane protein found only in the rough endoplasmic reticulum. The encoded protein is part of an N-oligosaccharyl transferase complex that links high mannose oligosaccharides to asparagine residues found in the Asn-X-Ser/Thr consensus motif of nascent polypeptide chains. This protein is similar in sequence to the yeast oligosaccharyl transferase subunit SWP1. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2008]