

## Product datasheet for **RC205872L2V**

### FAAP24 (NM\_152266) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	FAAP24 (NM_152266) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FAAP24
Synonyms:	C19orf40
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_152266
ORF Size:	645 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205872).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_152266.1</a>
RefSeq Size:	1130 bp
RefSeq ORF:	648 bp
Locus ID:	91442
UniProt ID:	<a href="#">Q9BTP7</a>
Cytogenetics:	19q13.11
MW:	23.7 kDa



[View online »](#)

**Gene Summary:**

FAAP24 is a component of the Fanconi anemia (FA) core complex (see MIM 227650), which plays a crucial role in DNA damage response (Ciccia et al., 2007 [PubMed 17289582]). [supplied by OMIM, Mar 2008]