

Product datasheet for **RC205877A1V**

Human IFI30 (NM_006332) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human IFI30 (NM_006332) AAV Particle
Tag:	Myc-DDK
Symbol:	IFI30
Synonyms:	GILT; IFI-30; IP-30; IP30
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC205877 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGACCCTGTCGCCACTTCTGCTGTTCTGCCACCGCTGCTGCTGCTGGACGTCCCCACGGCGGCGG
TGCAGGCGTCCCCTCTGCAAGCGTTAGACTTCTTTGGGAATGGGCCACCAGTTAACTACAAGACAGGCAA
TCTATACCTGCGGGGGCCCCGAAGAAGTCCAATGCACCGCTTGCAATGTGACCCTCTACTATGAAGCA
CTGTGCGGTGGCTGCCGAGCCTTCCTGATCCGGGAGCTTTCCCAACATGGCTGTTGGTCATGGAGATCC
TCAATGTCACGCTGGTGCCCTACGGAAACGCACAGGAACAAAATGTCAGTGGCAGGTGGGAGTTCAAGTG
CCAGCATGGAGAAGAGGAGTGCAAATTCACAAAGGTGGAGGCCCTGCGTGTGGATGAACCTGACATGGAG
CTAGCCTTCTGACCATTGTCTGCATGGAAGAGTTTGAGGACATGGAGAGAAGTCTGCCACTATGCCTGC
AGCTCTACGCCCCAGGGCTGTCGCCAGACACTATCATGGAGTGTGCAATGGGGACCGCGGCATGCAGCT
CATGCACGCCAACGCCAGCGGACAGATGCTCTCCAGCCACCACAGAGTATGTGCCCTGGGTCAACGTC
AATGGGAAACCTTGGAAGATCAGACCCAGCTCCTTACCCTTGCTGCCAGTTGTACCAGGGCAAGAAGC
CGGATGTCTGCCCTTCTCAACCAGCTCCCTCAGGAGTGTTCCTTCAAG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC205877 protein sequence Red=Cloning site Green=Tags(s) MTLSPLLLFLPPLLLLLDVPTAAVQASPLQALDFFGNGPPVNYKTGNLYLRGPKKSNAPLVNVTLYYEA LCGGCRAFLIRELFPTWLLVMEILNVTLPYGNAQEQNVSGRWEFKCQHGEEEKFNKVEACVLDELDME LAFLTIVCMEEFEDMERSLPLCLQLYAPGLSPDTIMECAMGDRGMQLMHANAQRDQLQPPHEYVPWTV NGKPLEDQTQLLTLVCQLYQGKKPDVCPSSSTSSLRVCFK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_006332
ORF Size:	750 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_006332.3</u>
RefSeq Size:	1034 bp
RefSeq ORF:	753 bp
Locus ID:	10437
UniProt ID:	<u>P13284</u>
Cytogenetics:	19p13.11
MW:	28 kDa