

Product datasheet for **RC205895L1V**

WDR11 (NM_018117) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	WDR11 (NM_018117) Human Tagged ORF Clone Lentiviral Particle
Symbol:	WDR11
Synonyms:	BRWD2; DR11; HH14; SRI1; WDR15
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_018117
ORF Size:	3672 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205895).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018117.10
RefSeq Size:	4746 bp
RefSeq ORF:	3675 bp
Locus ID:	55717
UniProt ID:	Q9BZH6
Cytogenetics:	10q26.12
Domains:	WD40
MW:	136.7 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the WD repeat protein family. WD repeats are minimally conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-asp (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. Members of this family are involved in a variety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene regulation. This gene is located in the chromosome 10q25-26 region, which is frequently deleted in gliomas and tumors of other tissues, and is disrupted by the t(10;19) translocation rearrangement in glioblastoma cells. The gene location suggests that it is a candidate gene for the tumor suppressor locus. [provided by RefSeq, Jul 2008]