

## Product datasheet for **RC205914L4V**

### **DHX15 (NM\_001358) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	DHX15 (NM_001358) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DHX15
Synonyms:	DBP1; DDX15; PRP43; PRPF43; PrPp43p
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001358
ORF Size:	2385 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205914).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001358.2</a>
RefSeq Size:	3019 bp
RefSeq ORF:	2388 bp
Locus ID:	1665
UniProt ID:	<a href="#">O43143</a>
Cytogenetics:	4p15.2
Domains:	DEAD, helicase_C, HA2
Protein Pathways:	Spliceosome



[View online »](#)

**MW:** 90.9 kDa

**Gene Summary:** The protein encoded by this gene is a putative ATP-dependent RNA helicase implicated in pre-mRNA splicing. [provided by RefSeq, Jul 2008]