

Product datasheet for **RC205928L3V**

hCG beta (CGB5) (NM_033043) Human Tagged ORF Clone Lentiviral Particle

Product data:


Product Type:	Lentiviral Particles
Product Name:	hCG beta (CGB5) (NM_033043) Human Tagged ORF Clone Lentiviral Particle
Symbol:	hCG beta
Synonyms:	CGB; HCG; hCGB
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_033043
ORF Size:	495 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205928).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_033043.1
RefSeq Size:	880 bp
RefSeq ORF:	498 bp
Locus ID:	93659
UniProt ID:	P01233
Cytogenetics:	19q13.33
Protein Families:	Druggable Genome
MW:	17.7 kDa



[View online »](#)

Gene Summary:

This gene is a member of the glycoprotein hormone beta chain family and encodes the beta 5 subunit of chorionic gonadotropin (CG). Glycoprotein hormones are heterodimers consisting of a common alpha subunit and an unique beta subunit which confers biological specificity. CG is produced by the trophoblastic cells of the placenta and stimulates the ovaries to synthesize the steroids that are essential for the maintenance of pregnancy. The beta subunit of CG is encoded by 6 genes which are arranged in tandem and inverted pairs on chromosome 19q13.3 and contiguous with the luteinizing hormone beta subunit gene. [provided by RefSeq, Jul 2008]

Product images:

[RC205928L3] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC205928L3V particle to overexpress human CGB5-Myc-DDK fusion protein.