

## Product datasheet for **RC205929L3V**

### **C11orf17 (AKIP1) (NM\_182901) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	C11orf17 (AKIP1) (NM_182901) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C11orf17
Synonyms:	BCA3; C11orf17
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_182901
ORF Size:	630 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205929).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_182901.2</a> , <a href="#">NP_878904.2</a>
RefSeq Size:	1528 bp
RefSeq ORF:	632 bp
Locus ID:	56672
Cytogenetics:	11p15.4
MW:	23.1 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a nuclear protein that interacts with protein kinase A catalytic subunit, and regulates the effect of the cAMP-dependent protein kinase signaling pathway on the NF-kappa-B activation cascade. Alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Oct 2011]