

## Product datasheet for **RC205932L4V**

### Myosin Light Chain 2 (MYL2) (NM\_000432) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Myosin Light Chain 2 (MYL2) (NM_000432) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Myosin Light Chain 2
Synonyms:	CMH10; MLC-2s/v; MLC2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_000432
ORF Size:	498 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205932).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_000432.1</a>
RefSeq Size:	855 bp
RefSeq ORF:	501 bp
Locus ID:	4633
UniProt ID:	<a href="#">P10916</a>
Cytogenetics:	12q24.11
Domains:	EFh



[View online »](#)

<b>Protein Pathways:</b>	Cardiac muscle contraction, Dilated cardiomyopathy, Focal adhesion, Hypertrophic cardiomyopathy (HCM), Leukocyte transendothelial migration, Regulation of actin cytoskeleton, Tight junction
<b>MW:</b>	18.8 kDa
<b>Gene Summary:</b>	Thus gene encodes the regulatory light chain associated with cardiac myosin beta (or slow) heavy chain. Ca <sup>+</sup> triggers the phosphorylation of regulatory light chain that in turn triggers contraction. Mutations in this gene are associated with mid-left ventricular chamber type hypertrophic cardiomyopathy. [provided by RefSeq, Jul 2008]