

Product datasheet for **RC205938A1V**

Human CRIP2 (NM_001312) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human CRIP2 (NM_001312) AAV Particle
Tag: Myc-DDK
Symbol: CRIP2
Synonyms: CRIP; CRP2; ESP1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC205938 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCCTCCAAATGCCCAAGTGCACAAAGACCGTGTACTTCGCCGAGAAGGTGAGCTCCCTGGGGAAGG
 ACTGGCACAAGTTCTGCCTCAAGTGCAGAGCGCTGCAGCAAGACGCTGACGCCCGGGGCCACGCCGAGCA
 TGACGGGAAGCCGTTCTGCCACAAGCCGTGCTACGCCACCCTGTTCCGACCCAAAGGCGTGAACATCGGG
 GGCGCGGGCTCCTACATCTACGAGAAGCCCTGGCGGAGGGGCCGAGGTACCGGCCCATCGAGGTCC
 CCGCGGCCGAGCAGAGGAGCGGAAGGCGAGCGGCCCCCGAAGGGGCCAGCAGAGCCTCCAGTGTAC
 CACTTTTACCGGGAGCCCAACACGTGCGCGCTGCAGCAAGAAGGTGTACTTCGCTGAGAAGGTGACG
 TCTCTGGGCAAGGATTGGCACCGGCCCTGCCTGCGCTGCGAGCGCTGCGGGAAGACACTGACCCCGGCG
 GGCACGCGGAGCACGACGCCAGCCCTACTGCCACAAGCCCTGCTATGGAATCCTCTTCGGACCCAAGGG
 AGTGAACACCGGTGCGGTGGGCAGCTACATCTATGACCGGGACCCCGAAGGCAAGGTCCAGCCC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC205938 protein sequence
 Red=Cloning site Green=Tags(s)

MASKCPKCDKTVYFAEKVSSLGKDWKFLKLCERCSKTLTPGGHAEHDGKPFCHKPCYATLFGPKGVNIG
 GAGSYIYEKPLAEGPQVTGPIEVPAARAEERKASGPPKGPSRASSVTTFTEPNTCPRCSKKVYFAEKVT
 SLGKDWHRPCLRCERCCKTLTPGGHAEHDGQPYCHKPCYGILFGPKGVNTGAVGSYIYDRDPEGKVQP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_001312
ORF Size:	624 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001312.2</u>
RefSeq Size:	1247 bp
RefSeq ORF:	627 bp
Locus ID:	1397
UniProt ID:	<u>P52943</u>
Cytogenetics:	14q32.33
MW:	22.5 kDa