

Product datasheet for **RC205938L4V**

CRIP2 (NM_001312) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CRIP2 (NM_001312) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CRIP2
Synonyms:	CRIP; CRP2; ESP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001312
ORF Size:	624 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205938).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001312.2
RefSeq Size:	1247 bp
RefSeq ORF:	627 bp
Locus ID:	1397
UniProt ID:	P52943
Cytogenetics:	14q32.33
Domains:	LIM
MW:	22.5 kDa



[View online »](#)

Gene Summary:

This gene encodes a putative transcription factor with two LIM zinc-binding domains. The encoded protein may participate in the differentiation of smooth muscle tissue. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2012]