

Product datasheet for **RC205962**

GFM2 (NM_032380) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GFM2 (NM_032380) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	GFM2
Synonyms:	EF-G2mt; EFG2; hEFG2; mEF-G 2; MRRF2; MST027; MSTP027; RRF; RRF2; RRF2mt
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC205962 ORF sequence

Red=Cloning site Blue=ORF Green=Tags(s)

 TTTTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

 ATGTTGACCAACTTGAGGATATTTGCAATGAGTCATCAGACAATACCCAGTGTGTATTAATAATATAT
GCTGCTATAAAATAAGAGCAAGTTTAAAAAGATTAAGCCACATGTGCCGCTTGGAAAGAAATTGCAGTTC
TCTACCAGGCTTAATAGGAAATGATATCAAATCCATTCTTCATCATCAATCCTCCCATAGCTAAATC
CGTAATATTGGAATTATGGCTCATATTGATGCAGGCAAACTACCACCACAGAAAGAATATTGTACTATT
CCGGATATACAAGATCACTGGGAGATGTTGATGATGGAGACACAGTGACAGATTCATGGCCCAAGAGCG
AGAAAGAGGCATTACTATTCAATCAGCTGCTGTTACATTTGATTGGAAAGGTTATAGAGTCAATCTAATT
GATACACCAGGTCATGTGGACTTTACCTTGGAGGTTGAACGGTGCCTAAGAGTGTGGATGGTGCAGTGG
CTGATTTTGATGCCTCTGCTGGTGTAGAGGCCAGACTCTCACAGTATGGAGGCAAGCTGATAAACACAA
TATACCTCGAATCTGTTTTTAAACAAGATGGACAAAAGTGGAGCAAGCTTTAAGTATGCAGTTGAAAGC
ATCAGAGAGAAGTTAAAGGCAAGCCTTTGCTTTTACAGTTACCAATTGGTGAAGCCAAAAGCTTTCAAAG
GAGTGGTGGATGTAGTAATGAAAGAAAACTTCTTTGGAATTGCAATTCAAATGATGGAAAAGACTTTGA
GAGAAAGCCCCTCTTGGAAATGAATGATCCTGAATTGCTGAAGGAAACAAGTGAAGCAAGGAATGCCTTA
ATTGAACAAGTTGCAGATTTGGATGATGAATTTGCTGACTTGGTTTTAGAAGATTTAGTGAGAATTTTG
ATTTGTTACCAGCTGAAAAGCTACAGACTGCAATACATAGAGTGACACTAGCTCAGACAGCAGTGCCTGT
GCTTTGTGGAAGTGCCCTGAAAAACAAAGGGATACAGCCCTTGTTAGATGCTGTTACTATGTACTTACCT
TCACCTGAAGAGCGTAACATGAATTTCTGCAGTGGTATAAGGATGACTTATGTGCATTGGCATTAAAG
TTCTCCATGACAAGCAGCGAGGACCACTGGTTTTTATGCGCATTTACTCAGGCACTATAAAACCCAGTT
GGCCATTCATAATATTAATGAAAAGTGCACGGAGAGAATAAGTCGCTCTGCTTTTGCCGTTTGCTGACCAA
CATGTAGAAATCCATTCTGACTGCTGGTAACATTGCTTTGACTGTTGGGCTTAAACATACTGCCACTG
GAGACACCATTGTCTCATCCAAGTCCAGTGCATTAGCTGCAGCTCGTAGAGCCGAACGGGAGGGAGAAAA
GAAGCACAGACAAAACAATGAAGCAGAGAGACTTTTATTGGCTGGAGTGGAGATTCCAGAACCTGTTTTT
TTCTGTACCATAGAACCCCATCACTGTCTAAGCAGCCAGATTTGGAACATGCGTTGAAATGTCTTCAGC
GTGAAGATCCCAGTTTGAAGTGAGGCTAGATCCTGACTCTGGACAACTGTTCTGTGTGGTATGGGGGA
GTTACATATAGAGATTATTCATGATCGAATCAAGAGGGAATATGGACTGGAGACCTATCTCGGGCCTCTC
CAGGTGGCATATCGAGAGACCATCCTAACTCAGTTTCGTGCCACAGATACCTTAGATAGAACTTTAGGAG
ACAAAAGGCATCTGTGACTGTAGAAGTGGGAGCAAGGCCAATTGAAACATCATCTGTTATGCCTGTGAT
TGAGTATGCTGAAAGTATCAATGAAGGCCTTTTGAAGGCTCACAAGAGGCCATTGAAAAATGGAATTCAC
AGCGCATGTCTCCAAGGACCATTTGCTTGGATCCCCAATTCAGGATGTAGCAATTACTTTACATTCCTCTGA
CAATTCATCCTGGCACCTCCACAACTATGATTTCTGCCTGTGTCTCAAGATGCGTGCAAAAGGCTCTGAA
GAAAGCTGATAAGCAAGTTTGGAGCCTCTGATGAATCTTGAGGTTACAGTAGCTAGAGATTATCTCAGC
CCTGTCTGGCAGATCTGGCACAAGAAGAGGAAACATTCAGGAAATTCAGACTCGCCAGGACAACAAAG
TTGTTATTGGATTGTTCCTTAGCAGGAATTATGGGTTATTCACTGTGCTTCGAACGCTAACATCAGG
CTCAGCTACTTTTGCTTAGAACTATCTACTTATCAAGCCATGAATCCTCAAGATCAAAATACACTGCTC
AACCGGAGAAGTGGTTTGACC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC205962 protein sequence
 Red=Cloning site Green=Tags(s)


MLTNLRIFAMSHQTIPSVYINNICCCKIRASLRLKPHVPLGRNCSSLPGLIGNDIKSIHSIINPPIAKI
 RNIGIMAHIDAGKTTTTTERILYYSGYTRSLGDVDDGDTVDFMAQERERGITIQSAAVTFDWKGYRVNLI
 DTPGHVDFTLVERCLRVLDGAVAVFDASAGVEAQTLTVWRQADKHNIPRICFLNKMDKTGASFYAVES
 IREKLKAKPLLLQLPIGEAKTFKGVVDVVMKEKLLWNCNSNDGKDFERKPLEMNDPELLKETTEARNAL
 IEQVADLDDEFADLVLEEFSENFDLLPAEKLQTAIHRVTLAQTAVPVLCGSALKNKGIQPLLDVATMYLP
 SPEERNYEFLQWYKDDLALAFKVLHDKQRGPLVFMRIYSGTIKPQLAIHNINGNCTERISRLLLPFADQ
 HVEIHSLTAGNIALTVGLKHTATGDTIVSSKSSALAAARRAEREKKEKHRQNEAERLLLAGVEIPEPVF
 FCTIEPPSLSKQPDLEHALKCLQREDPSLKVRDPDSGQTVLCGMGELHIEIIHDRIKREYGLETYLGPL
 QVAYREILNSVRATDTLDRTLGDKRHLVTVEVGARPIETSSVMPVIEYAESINEGLLKVSQEAIENGIH
 SACLQGPLLGSPIQDVAITLHSLTIHPGTSTTMISACVSRVQKALKKADKQVLEPLMNLEVTVARDYLS
 PVLADLAQRNGNIQEIQTRQDNKVVGIVPLAGIMGYSTVLRTLTSGSATFALELSTYQAMNPQDQNTLL
 NRRSGLT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6695_b12.zip

Restriction Sites: SgfI-MluI


Cloning Scheme:




ACCN:	NM_032380
ORF Size:	2331 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq Size:	3264 bp
RefSeq ORF:	2340 bp
Locus ID:	84340
UniProt ID:	<u>Q969S9</u>
Cytogenetics:	5q13.3
Domains:	EFG_C, GTP_EFTU, GTP_EFTU_D2, EFG_IV
MW:	86.2 kDa
Gene Summary:	Eukaryotes contain two protein translational systems, one in the cytoplasm and one in the mitochondria. Mitochondrial translation is crucial for maintaining mitochondrial function and mutations in this system lead to a breakdown in the respiratory chain-oxidative phosphorylation system and to impaired maintenance of mitochondrial DNA. This gene encodes one of the mitochondrial translation elongation factors, which is a GTPase that plays a role at the termination of mitochondrial translation by mediating the disassembly of ribosomes from messenger RNA . Its role in the regulation of normal mitochondrial function and in disease states attributed to mitochondrial dysfunction is not known. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2013]


Product images:



Circular map for RC205962



Western blot validation of overexpression lysate (Cat# [LY410173]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC205962 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified GFM2 protein (Cat# [TP305962]). The protein was produced from HEK293T cells transfected with GFM2 cDNA clone (Cat# RC205962) using MegaTran 2.0 (Cat# [TT210002]).