

Product datasheet for **RC205970L4V**

TLL2 (NM_031949) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TLL2 (NM_031949) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TLL2
Synonyms:	C6orf104; dj366N23.3; NYD-TSPG
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031949
ORF Size:	1776 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205970).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031949.4
RefSeq Size:	2877 bp
RefSeq ORF:	1779 bp
Locus ID:	83887
UniProt ID:	Q9BWV7
Cytogenetics:	6q27
MW:	67.3 kDa



[View online »](#)

Gene Summary:

Probable tubulin polyglutamylase that forms polyglutamate side chains on tubulin. Probably acts when complexed with other proteins (By similarity).[UniProtKB/Swiss-Prot Function]