

Product datasheet for **RC206001L3V**

SERINC1 (NM_020755) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SERINC1 (NM_020755) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SERINC1
Synonyms:	TDE1L; TDE2; TMS-2; TMS2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_020755
ORF Size:	1359 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206001).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_020755.2
RefSeq Size:	3144 bp
RefSeq ORF:	1362 bp
Locus ID:	57515
UniProt ID:	Q9NRX5
Cytogenetics:	6q22.31
Domains:	TMS_TDE
Protein Families:	Transmembrane



[View online »](#)

MW: 50.5 kDa

Gene Summary: Enhances the incorporation of serine into phosphatidylserine and sphingolipids.
[UniProtKB/Swiss-Prot Function]