

## Product datasheet for RC206004L2V

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## CRM1 (XPO1) (NM\_003400) Human Tagged ORF Clone Lentiviral Particle

**Product data:** 

Product Type: Lentiviral Particles

Product Name: CRM1 (XPO1) (NM\_003400) Human Tagged ORF Clone Lentiviral Particle

Symbol: XPO1

Synonyms: CRM-1; CRM1; emb; exp1

Mammalian Cell

Selection:

None

**Vector:** pLenti-C-mGFP (PS100071)

Tag: mGFP

**ACCN:** NM\_003400 **ORF Size:** 3213 bp

**ORF Nucleotide** 

The ORF insert of this clone is exactly the same as(RC206004).

Sequence:

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through

naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**RefSeg:** NM 003400.3

 RefSeq Size:
 4830 bp

 RefSeq ORF:
 3216 bp

 Locus ID:
 7514

 UniProt ID:
 014980

 Cytogenetics:
 2p15

Domains: IBN NT

**Protein Families:** Druggable Genome





**MW:** 123.8 kDa

**Gene Summary:** This cell-cycle-regulated gene encodes a protein that mediates leucine-rich nuclear export

signal (NES)-dependent protein transport. The protein specifically inhibits the nuclear export of Rev and U snRNAs. It is involved in the control of several cellular processes by controlling the localization of cyclin B, MPAK, and MAPKAP kinase 2. This protein also regulates NFAT and

AP-1. [provided by RefSeq, Jan 2015]