

Product datasheet for **RC206022A1V**

Human Kallikrein 11 (KLK11) (NM_006853) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Kallikrein 11 (KLK11) (NM_006853) AAV Particle
Tag:	Myc-DDK
Symbol:	Kallikrein 11
Synonyms:	PRSS20; TLSP
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC206022 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAGGATTCTGCAGTTAATCCTGCTTCTGGCAACAGGGCTTGTAGGGGAGAGACCAGGATCATCA
AGGGGTTTCAGTGAAGCCTCACTCCAGCCCTGGCAGGCAGCCCTGTTGAGAAGACGCGGCTACTCTG
TGGGGCAGCGCTCATCGCCCCAGATGGCTCCTGACAGCAGCCCACTGCCTCAAGCCCCGCTACATAGTT
CACCTGGGGCAGCACAACCTCCAGAAGGAGGAGGGCTGTGAGCAGACCCGGACAGCCACTGAGTCCTTCC
CCCACCCCGGCTTCAACAACAGCCTCCCAACAAGACCACCGCAATGACATCATGCTGGTGAAGATGGC
ATCGCCAGTCTCCATCACCTGGGCTGTGCGACCCCTCACCTCTCCTCACGCTGTGTCACTGCTGGCACC
AGCTGCCTCATTTCCGGCTGGGGCAGCACGTCCAGCCCCAGTTACGCTGCCTCACACCTTGCGATGCG
CCAACATCACCATCATTGAGCACCAGAAGTGTGAGAAGCCTACCCCGGCAACATCACAGACACCATGGT
GTGTGCCAGCGTGCAGGAAGGGGGCAAGGACTCCTGCCAGGGTGAATCCGGGGGCCCTCTGGTCTGTAAC
CAGTCTCTTCAAGGCATTATCTCTGGGGCCAGGATCCGTGTGCGATCACCCGAAAGCCTGGTGTCTACA
CGAAAGCTGCAAAATATGTGGACTGGATCCAGGAGACGATGAAGAACAAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC206022 protein sequence Red=Cloning site Green=Tags(s) MRILQLILLALATGLVGGETRIIKGFECKPHSQPWQAALFEKTRLLCGATLIAPRWLLTAAHCLKPRYIV HLGQHNLQKEEGCEQTRTATESFPHPGFNNSLPNKDHRNDIMLVKMASPVSTITWAVRPLTLSSRCVTAGT SCLISGWGSTSSPQLRLPHTLRANITIIIEHQKCENAYPGNITDTMVCASVQEGGKDSCQGDSSGGLVCN QSLQGIISWGQDPCAITRKPGVYTKVCKYVDWIQETMKNN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_006853
ORF Size:	750 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_006853.2</u>
RefSeq Size:	1207 bp
RefSeq ORF:	753 bp
Locus ID:	11012
UniProt ID:	<u>Q9UBX7</u>
Cytogenetics:	19q13.41
MW:	27.5 kDa