

Product datasheet for **RC206027L3V**

VPS45A (VPS45) (NM_007259) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	VPS45A (VPS45) (NM_007259) Human Tagged ORF Clone Lentiviral Particle
Symbol:	VPS45A
Synonyms:	H1; H1VPS45; SCN5; VPS45A; VPS45B; VPS54A; VSP45; VSP45A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_007259
ORF Size:	1710 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206027).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007259.3
RefSeq Size:	2824 bp
RefSeq ORF:	1713 bp
Locus ID:	11311
UniProt ID:	Q9NRW7
Cytogenetics:	1q21.2
Domains:	Sec1
Protein Pathways:	Endocytosis



[View online »](#)

MW: 65.1 kDa

Gene Summary: Vesicle mediated protein sorting plays an important role in segregation of intracellular molecules into distinct organelles. Genetic studies in yeast have identified more than 40 vacuolar protein sorting (VPS) genes involved in vesicle transport to vacuoles. This gene is a member of the Sec1 domain family, and shows a high degree of sequence similarity to mouse, rat and yeast Vps45. The exact function of this gene is not known, but its high expression in peripheral blood mononuclear cells suggests a role in trafficking proteins, including inflammatory mediators. Multiple alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jul 2013]