

Product datasheet for **RC206028L4V**

KRR1 (NM_007043) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	KRR1 (NM_007043) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KRR1
Synonyms:	HRB2; RIP-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007043
ORF Size:	1143 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206028).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007043.6
RefSeq Size:	3357 bp
RefSeq ORF:	1146 bp
Locus ID:	11103
UniProt ID:	Q13601
Cytogenetics:	12q21.2
Domains:	KH
MW:	43.5 kDa



[View online »](#)

Gene Summary: Required for 40S ribosome biogenesis. Involved in nucleolar processing of pre-18S ribosomal RNA and ribosome assembly (By similarity).[UniProtKB/Swiss-Prot Function]