

Product datasheet for **RC206037L4V**

ST5 (DENND2B) (NM_213618) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ST5 (DENND2B) (NM_213618) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DENND2B
Synonyms:	HTS1; p126; ST5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_213618
ORF Size:	3411 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206037).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_213618.1 , NP_998783.1
RefSeq Size:	4371 bp
RefSeq ORF:	3414 bp
Locus ID:	6764
UniProt ID:	P78524
Cytogenetics:	11p15.4
MW:	126.4 kDa



[View online »](#)

Gene Summary:

This gene was identified by its ability to suppress the tumorigenicity of HeLa cells in nude mice. The protein encoded by this gene contains a C-terminal region that shares similarity with the Rab 3 family of small GTP binding proteins. This protein preferentially binds to the SH3 domain of c-Abl kinase, and acts as a regulator of MAPK1/ERK2 kinase, which may contribute to its ability to reduce the tumorigenic phenotype in cells. Three alternatively spliced transcript variants of this gene encoding distinct isoforms are identified. [provided by RefSeq, Jul 2008]