

Product datasheet for **RC206046A1V**

Human COX5A (NM_004255) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human COX5A (NM_004255) AAV Particle
Tag: Myc-DDK
Symbol: COX5A
Synonyms: COX; COX-VA; MC4DN20; VA
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC206046 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCTGGGCGCCGCTCTCCGCCGCTGCGCTGTGGCCGCAACCCGGGCGCACCTCGAGGCCTCCTGC
ACTCCGCCCGGACCCCGGCCCGCGGTGGCTATCCAGTCAGTTCGCTGCTATCCCATGGGTACAGGA
GACAGATGAGGAGTTTGATGCTCGCTGGGTAACATACTTCAACAAGCCAGATATAGATGCCTGGGAATTG
CGTAAAGGGATAAACACACTTGTTACCTATGATATGGTCCAGAGCCCCAAATCATTGATGCTGCTTTGC
GGGCATGCAGACGGTTAAATGATTTTGCTAGTACAGTTCGTATCCTAGAGTTGTTAAGGACAAAGCAGG
ACCTCATAAGGAAATCTACCCCTATGTCATCCAGGAACCTAGACCAACTTTAAATGAAGTGGGAATCTCC
ACTCCGAGGAACTGGGCCTTGACAAAGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC206046 protein sequence
Red=Cloning site Green=Tags(s)

MLGAALRRCAVAATTRADPRGLLHSARTPGPAVAIQSVRCYSHGSQETDEEFDARWVTFYFNKPDIDAWEL
RKGINTLVTYDMVPEPKIIDAALRACRRLNDFASTVRILEVVKDKAGPHKEIYPYVIQELRPTLNELGIS
TPEELGLDKV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online »](#)

ACCN:	NM_004255
ORF Size:	450 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_004255.2</u>
RefSeq Size:	784 bp
RefSeq ORF:	453 bp
Locus ID:	9377
UniProt ID:	<u>P20674</u>
Cytogenetics:	15q24.2
MW:	16.8 kDa