

Product datasheet for **RC206052L4V**

PRPS2 (NM_002765) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PRPS2 (NM_002765) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PRPS2
Synonyms:	PRSII
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_002765
ORF Size:	954 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206052).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002765.3
RefSeq Size:	2509 bp
RefSeq ORF:	957 bp
Locus ID:	5634
UniProt ID:	P11908
Cytogenetics:	Xp22.2
Domains:	Pribosyltran
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Metabolic pathways, Pentose phosphate pathway, Purine metabolism

MW: 34.6 kDa

Gene Summary: This gene encodes a phosphoribosyl pyrophosphate synthetase that plays a central role in the synthesis of purines and pyrimidines. The encoded protein catalyzes the synthesis of 5-phosphoribosyl 1-pyrophosphate from ATP and D-ribose 5-phosphate. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Mar 2010]