

## Product datasheet for **RC206063L2V**

### **GUCY1A3 (GUCY1A1) (NM\_000856) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	GUCY1A3 (GUCY1A1) (NM_000856) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GUCY1A1
Synonyms:	GC-S-alpha-1; GC-SA3; GCS-alpha-3; GUC1A3; GUCA3; GUCSA3; GUCY1A3; MYMY6
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_000856
ORF Size:	2070 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206063).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_000856.2</a>
RefSeq Size:	9454 bp
RefSeq ORF:	2073 bp
Locus ID:	2982
UniProt ID:	<a href="#">Q02108</a>
Cytogenetics:	4q32.1
Domains:	CYcC
Protein Families:	Druggable Genome



[View online »](#)

**Protein Pathways:** Gap junction, Long-term depression, Purine metabolism, Vascular smooth muscle contraction

**MW:** 77.5 kDa

**Gene Summary:** Soluble guanylate cyclases are heterodimeric proteins that catalyze the conversion of GTP to 3',5'-cyclic GMP and pyrophosphate. The protein encoded by this gene is an alpha subunit of this complex and it interacts with a beta subunit to form the guanylate cyclase enzyme, which is activated by nitric oxide. Several transcript variants encoding a few different isoforms have been found for this gene. [provided by RefSeq, Jan 2012]