

Product datasheet for **RC206069L1V**

Aquaporin 5 (AQP5) (NM_001651) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Aquaporin 5 (AQP5) (NM_001651) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Aquaporin 5
Synonyms:	AQP-5; PPKB
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001651
ORF Size:	795 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206069).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001651.1
RefSeq Size:	1840 bp
RefSeq ORF:	798 bp
Locus ID:	362
UniProt ID:	P55064
Cytogenetics:	12q13.12
Protein Families:	Druggable Genome, Transmembrane
MW:	28.3 kDa



[View online »](#)

Gene Summary:

Aquaporin 5 (AQP5) is a water channel protein. Aquaporins are a family of small integral membrane proteins related to the major intrinsic protein (MIP or AQP0). Aquaporin 5 plays a role in the generation of saliva, tears and pulmonary secretions. AQP0, AQP2, AQP5, and AQP6 are closely related and all map to 12q13. [provided by RefSeq, Jul 2008]