

Product datasheet for **RC206099L3V**

GTSF1 (NM_144594) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GTSF1 (NM_144594) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GTSF1
Synonyms:	Cue110; FAM112B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_144594
ORF Size:	501 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206099).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_144594.1 , NP_653195.1
RefSeq Size:	822 bp
RefSeq ORF:	504 bp
Locus ID:	121355
UniProt ID:	Q8WW33
Cytogenetics:	12q13.13
MW:	19.2 kDa



[View online »](#)

Gene Summary: Required for spermatogenesis and is involved in the suppression of retrotransposon transcription in male germ cells.[UniProtKB/Swiss-Prot Function]