

Product datasheet for **RC206107L4V**

TTC18 (CFAP70) (NM_145170) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TTC18 (CFAP70) (NM_145170) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TTC18
Synonyms:	TTC18
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145170
ORF Size:	3363 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206107).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145170.2 , NP_660153.2
RefSeq Size:	3714 bp
RefSeq ORF:	3366 bp
Locus ID:	118491
Cytogenetics:	10q22.2
MW:	125.7 kDa
Gene Summary:	Axoneme-binding protein that plays a role in the regulation of ciliary motility and cilium length.[UniProtKB/Swiss-Prot Function]



[View online »](#)