

Product datasheet for **RC206114A1V**

Human ASF1B (NM_018154) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human ASF1B (NM_018154) AAV Particle
Tag: Myc-DDK
Symbol: ASF1B
Synonyms: CIA-II
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC206114 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCAAGGTGTCGGTGCTGAACGTGGCGGTCTGGAGAACCCGAGCCCTTTCCACAGCCCTTCCGGT
TCGAGATCAGCTTCGAGTGCAGTGAAGCCCTGGCGGACGACCTGGAGTGAAGATCATTTATGTTGGCTC
GGCTGAGAGTGAGGAATTTGATCAGATCCTAGACTCGGTGCTGGTGGGCCCTGTGCCAGCAGGGAGACAC
ATGTTTGTCTTTCAGGCCGACGCCCCAACCCATCCCTCATCCAGAGACTGATGCCGTGGGTGTGACTG
TGGTCCTCATCACCTGCACCTACCATGGACAGGAGTTTCATCCGAGTGGGCTACTACGTCAACAACGAGTA
CCTCAACCCCTGAGCTGCGTGAGAACCCGCCCATGAAGCCAGATTCTCCAGCTCCAGCGGAACATCTTG
GCCTCGAACCCCGGGTGACCCGCTTCCATATCAACTGGGACAACAACATGGACAGGCTGGAGGCCATAG
AGACCCAGGACCCCTCCCTGGGCTGCGGCCTCCCACTCAACTGCACTCCTATCAAGGGCTTGGGGCTCCC
TGGCTGCATCCCTGGCTCCTCCCTGAGAACTCCATGGACTGCATC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC206114 protein sequence
Red=Cloning site Green=Tags(s)

MAKVSVLNVAVLENPSPFHSPFRFEISFECSEALADDLEWKIIYVGSASEEFDQILDSVLVGPVPAGRH
MFVFQADAPNPSLIPETDAVGTVVLITCTYHGQEFIRVGYYVNNEYLNPELRENPPMKPDFSQLQRNII
ASNPRVTRFHINWDNMDRLEAIETQDPSLGCGLPLNCTPIKGLGLPGCIPGLLPENSMDCI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_018154
ORF Size:	606 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_018154.2</u>
RefSeq Size:	1746 bp
RefSeq ORF:	609 bp
Locus ID:	55723
UniProt ID:	<u>Q9NVP2</u>
Cytogenetics:	19p13.12
MW:	22.4 kDa