

## Product datasheet for **RC206118L4V**

### MFSD8 (NM\_152778) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	MFSD8 (NM_152778) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MFSD8
Synonyms:	CCMD; CLN7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_152778
ORF Size:	1554 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206118).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_152778.1</a>
RefSeq Size:	4562 bp
RefSeq ORF:	1557 bp
Locus ID:	256471
UniProt ID:	<a href="#">Q8NHS3</a>
Cytogenetics:	4q28.2
Domains:	sugar_tr
Protein Families:	Transmembrane



[View online »](#)

**Protein Pathways:** Lysosome

**MW:** 57.6 kDa

**Gene Summary:** This gene encodes a ubiquitous integral membrane protein that contains a transporter domain and a major facilitator superfamily (MFS) domain. Other members of the major facilitator superfamily transport small solutes through chemiosmotic ion gradients. The substrate transported by this protein is unknown. The protein likely localizes to lysosomal membranes. Mutations in this gene are correlated with a variant form of late infantile-onset neuronal ceroid lipofuscinoses (vLINCL). [provided by RefSeq, Oct 2008]