

Product datasheet for **RC206125A1V**

Human SPSB3 (NM_080861) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human SPSB3 (NM_080861) AAV Particle
Tag:	Myc-DDK
Symbol:	SPSB3
Synonyms:	C16orf31; SSB3
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC206125 representing NM_080861 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCCAGACGCCCCCGGAACAGCAGGGCCTGGCACTTCGTCTGAGTGACGCCCGGAGACGCAGATG
CCCCGGCCGTGGCTCTAGCAGGCTCCACTAACTGGGGCTACGACTCTGATGGGCAGCACAGCGACTCGGA
CTCCGACCCCGAGTACTCCACGCTGCCGCCATCCATCCCAAGTGGGTGCCGTGACCGGCGAGTCTTTC
TGTGACTGTGCTGGGCAGAGCGAGGCCCTCTTCTGTAGCAGCCTGCACTCGGCCACCGGGCAGGGACT
GCCGCTGCGGAGAGGAAGACGAGTATTTGCACTGGGTCTGGGATGACTTAAATAAGTCATCAGCCACCT
GCTGAGCTGTGACAACCGTAAGGTCAGCTTCCACATGGAGTACAGCTGCGGCACAGCGGCCATCCGGGGC
ACCAAGGAGCTGGGGAGGGCCAGCACTTCTGGGAGATCAAGATGACCTCTCCCGTCTACGGCACCGACA
TGATGGTGGGCATCGGGACGTCGGATGTGGACCTGGACAAATACCGCCACACGTTCTGCAGCCTGCTGGG
CAGGGATGAGGACAGCTGGGGCCTCTCTACACGGGCTCCTCCACCACAAGGGCGACAAGACCAGCTTC
TCGTCGCGGTTCCGGCAGGGCTCCATCATTGGCGTGCACCTGGACACCTGGCACGGCACACTCACCTTTT
TCAAGAACAGGAAGTGTATAGGTGTGGCAGCCACCAAGCTGCAGAACAGAGATTCTACCCGATGGTGTG
CTCCACGGCGGCGGAGCAGCATGAAGGTACCCGCTCCTGTGCCAGCGCCACTTCCCTCCAGTACCTG
TGCTGCCACCGCTGCGCCAGCTGCGGCAGACTCGGGAGACACGCTGGAGGGTCTGCCGCTGCCCGCGG
GCCTCAAGCAGGTGCTACACAACAAGCTGGGCTGGGTCTGAGCATGAGTTGCAGCCGCCGAAGGCTCC
AGTGTCCGATCCCGAGGAGCGACCTCCGCCACCACAGCAGTCGCGAGCCTCGGCCCTGCCAGAGGAAG
CGCTGCCCGCGGACC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC206125 representing NM_080861 Red=Cloning site Green=Tags(s) MARRPRNSRAWHFVLSAARRDADARAVALAGSTNWDGSDGQHSDDSDPEYSTLPSPISAVPVTGESF CDCAGQSEASFCSLHSAHRGRDCRCGEEDFYFDWVDDLKSSATLLSCDNRKVSFHMESCGTAAIRG TKELGEGQHFWEIKMTSPVYGDMVMGIGTSDVDLDKYRHTFCSLGRDEDSWGLSYTGLLHHKGDKTSF SSRFGQGSIIIGVHLDTWGTLTFFKNRKCIGVAATKLQNKRFYPMVCSTAARSSMKVTRSCASATSLQYL CCHRLRQLRPDSGDTLEGLPLPPGLKQVLHNKLGWVLSMSCSRRKAPVSDPQAATSAHSSREPRPCQRK RCRRT TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_080861
ORF Size:	1065 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_080861.3
RefSeq Size:	1553 bp
RefSeq ORF:	1068 bp
Locus ID:	90864
UniProt ID:	Q6PJ21
Cytogenetics:	16p13.3
MW:	39.2 kDa