

Product datasheet for **RC206125L4V**

SPSB3 (NM_080861) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SPSB3 (NM_080861) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SPSB3
Synonyms:	C16orf31; SSB3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_080861
ORF Size:	1065 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206125).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_080861.3
RefSeq Size:	1553 bp
RefSeq ORF:	1068 bp
Locus ID:	90864
UniProt ID:	Q6PJ21
Cytogenetics:	16p13.3
Domains:	SPRY
Protein Families:	Druggable Genome


[View online »](#)

MW: 39.2 kDa

Gene Summary: May be a substrate recognition component of a SCF-like ECS (Elongin BC-CUL2/5-SOCS-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.[UniProtKB/Swiss-Prot Function]