

Product datasheet for **RC206129A1V**

Human SUNC1 (SUN3) (NM_152782) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human SUNC1 (SUN3) (NM_152782) AAV Particle
Tag:	Myc-DDK
Symbol:	SUNC1
Synonyms:	SUNC1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC206129 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAGTGGAAAAACAAAGGCAAGAAGGGCTGCCATGTTTTTTAGACGTTGCTCTGAAGACGCCAGCGGTA
GCGCCAGTGGCAATGCTTTGTTATCAGAGGACGAAAATCCTGATGCGAATGGGGTAACTCGATCATGGAA
GATTATTCTAAGTACAATGCTTACACTGACTTTTCTTCTGTAGGACTCCTAAATCATCAGTGGCTTAAA
GAAACAGATGTTCTCAGAAATCCAGACAATTATATGCCATAATTGCAGAATATGGTTCAAGGCTTTATA
AATATCAGGCCAGACTTCGTATGCCTAAAGAGCAACTGGAACCTTTAAAGAAGGAAAGCCAGAATCTGGA
AAACAATTTTCGTCAAATTCATTTTTTGGTCGAACAAATAGATGTCCTGAAGGCATTGCTAAGAGATATG
AAGGATGGTATGGACAATAATCACAACCTGGAACACCCATGGAGACCCTGTGGAGGACCCGACCACACAG
AGGAAGTGTCAAACCTGGTCAATTATGTACTTAAAAAGTTGAGAGAAGACCAAGTCGAGATGGCTGATTA
TGCCCTGAAGTCGGCCGAGCCTCCATCATTGAAGCTGGGACCTCAGAAAGTTATAAAAAATAATAAGCA
AAATTGTACTGGCATGGGATAGGTTTCCATAATCATGAAATGCCTCCAGATATTATTCTCAGCCGGATG
TCTACCCTGGAAAGTGTGGGCTTTCCAGGTTCCAGGGTCATACCCTAATCAAGCTTGCTACAAAGAT
CATACCAACTGCTGTTACCATGGAGCACATCTCAGAGAAGGTGTCTCCGTCAGGAAACATCTCCAGTGCA
CCCAAGGAATTTTCTGTCTATGGCATCACAAAAAATGTGAAGGAGAAGAAATTTTCTAGGTCAAGTTTA
TATATAACAAAACAGGAACCACCGTTCAAACATTTGAAGTCCAGCATGCAGTTTCTGAATATTTATTATG
TGTGAAACTTAATATCTTTAGCAACTGGGGACACCCGAAGTATACTTGTTTATATCGATTCAAGGTCCAT
GGCACACCAGGCAAGCACATC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence:	>RC206129 protein sequence Red=Cloning site Green=Tags(s) MSGKTKARRAAMFFRRCSEDASGSASGNALLSEDPDANGVTRSWKIILSTMLTLFLLVGLLNHQWLK ETDVPQKSRQLYAIIEYGSRLYKYQARLRMPKEQLELLKKESQNLNNFRQILFLVEQIDVLKALLRDM KDGMDNNHNWNTHGDPVEDPDHTEEVSNLVNYVLKKLREDQVEMADYALKSAGASIEAGTSESYKNNKA KLYWHGIGFLNHEMPDIIILQPDVYPGKCWAFPGSQGHTLIKLATKIIPAVTMEHISEKVSPSGNISSA PKEFSVYGITKKCEGEEIFLGQFIYNKTGTTVQTFELQHAVSEYLLCVKLNIFSNWGHPKYTCLYRFRVH GTPGKHI TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_152782
ORF Size:	1071 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_152782.3</u> , <u>NP_689995.3</u>
RefSeq Size:	1368 bp
RefSeq ORF:	1074 bp
Locus ID:	256979
UniProt ID:	<u>Q8TAQ9</u>
Cytogenetics:	7p12.3
MW:	40.5 kDa