

## Product datasheet for **RC206130L3V**

### **DIXDC1 (NM\_033425) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	DIXDC1 (NM_033425) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DIXDC1
Synonyms:	CCD1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_033425
ORF Size:	1416 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206130).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_033425.3</a> , <a href="#">NP_219493.1</a>
RefSeq Size:	5496 bp
RefSeq ORF:	1419 bp
Locus ID:	85458
UniProt ID:	<a href="#">Q155Q3</a>
Cytogenetics:	11q23.1
MW:	53.9 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a positive regulator of the Wnt signaling pathway. The encoded protein is found associated with gamma tubulin at the centrosome. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jun 2013]