

Product datasheet for **RC206147L4V**

FAM71A (NM_153606) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FAM71A (NM_153606) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FAM71A
Synonyms:	GARI-L4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_153606
ORF Size:	1782 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206147).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_153606.1
RefSeq Size:	2332 bp
RefSeq ORF:	1785 bp
Locus ID:	149647
UniProt ID:	Q8IYT1
Cytogenetics:	1q32.3
MW:	63.3 kDa


[View online »](#)

Gene Summary:

This gene encodes a protein that interacts with the Rab2B small GTPase and may be important for integrity of the Golgi body. A knockdown of this gene induces fragmentation of the Golgi, similar to the effect seen with a knockdown of the Rab2B small GTPase. The encoded protein has an N-terminal Rab-binding domain specific for Rab2B. [provided by RefSeq, Feb 2017]