

Product datasheet for **RC206164L4V**

DNER (NM_139072) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DNER (NM_139072) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DNER
Synonyms:	bet; UNQ26
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_139072
ORF Size:	2211 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206164).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_139072.3
RefSeq Size:	3272 bp
RefSeq ORF:	2214 bp
Locus ID:	92737
UniProt ID:	Q8NFT8
Cytogenetics:	2q36.3
Domains:	EGF_CA, EGF, EGF
Protein Families:	Druggable Genome, Transmembrane


[View online »](#)

MW: 78.9 kDa

Gene Summary: Activator of the NOTCH1 pathway. May mediate neuron-glia interaction during astrocytogenesis (By similarity).[UniProtKB/Swiss-Prot Function]