

Product datasheet for **RC206173L2V**

RAB11FIP5 (NM_015470) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RAB11FIP5 (NM_015470) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RAB11FIP5
Synonyms:	gaf-1; GAF1; pp75; rab11-FIP5; RIP11
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_015470
ORF Size:	1956 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC206173).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015470.1
RefSeq Size:	4340 bp
RefSeq ORF:	1962 bp
Locus ID:	26056
UniProt ID:	Q9BXF6
Cytogenetics:	2p13.2
Protein Pathways:	Endocytosis
MW:	70.2 kDa



[View online »](#)

Gene Summary:

Rab effector involved in protein trafficking from apical recycling endosomes to the apical plasma membrane. Involved in insulin granule exocytosis. May regulate V-ATPase intracellular transport in response to extracellular acidosis.[UniProtKB/Swiss-Prot Function]